CRYOPLUS-2.0
CRYOPLUS-2.1

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CRYOPLUS-2.0
CRYOPLUS-2.1

These gloves are made with a special composite fabric bonded to a hydrophilic membrane. The interior features another full, breathable membrane, which does not allow liquids to penetrate inside. Protection from contact cold is guaranteed by an internal lining in soft polyester fleece.

CE CERTIFICATION

CE certification issued by: Ricotest Srl – Notified Lab No. 0498
Subject to specific liquid nitrogen resistance test (-195.82°Cc), by immersion in contact for 1 minute and subsequent RCT test.

N.B.: the performance level concerns the palm side, including fingers; nevertheless, the uniformity of materials and processing makes the back of the glove also protective. The following marking can be found on the label sewn into the inside of the gloves:

	CE marking, testifying to the conformity of the gloves to the essential health and safety requirements of Regulation (EU) 2016/425
0302	Identification no. of the Notified Body performing the inspection of the PPE manufactured pursuant (MODEL C2)
CRYOKIT	Manufacturer's registered trade mark
CRYOPLUS-2.0/2.1	Glove model
10 (example)	Size
	Pictogram referring to the information note.
abcde	EN388+pictogram relating to the mechanical risks explained below: 2(a) 3(b) 2(c) 4(d) x(e) This pictogram indicates that the gloves were designed for mechanical risks and the numbers at the foot indicate the levels obtained during laboratory tests (see explanation on the last page).
abc	EN 511 + related pictogram for protection against cold explained below: (a) 2 = level obtained against convective cold (0.193 m ² K/W) (b) 2 = level obtained against contact cold (0.1563 m ² K/W) (c) 1 = level obtained against water penetration (no penetration)

DESCRIPTION and CHARACTERISTICS

These gloves are category III Personal Protective Equipment (PPE) and are therefore subject to annual. They bear the CE marking since they conform to the requirements set by Regulation (EU) 2016/425 and have been certified by Notified Body. These gloves also comply with the requirements of harmonised technical standards EN 420:2003+A1:2009, EN 388:2016+A1 2018,, EN 511:2006 and EN ER07 (2004) only CRYOKIT HP and HPS.

RECOMMENDED USES/APPLICATIONS

Gloves intended for protection against the risks deriving from liquefied gases (cold burns and frostbite from the intense cold, leaking liquid, splashes) during loading, storage and decanting. Loading/storage/decanting of cryogenic liquefied gases by staff in charge of handling, contact with cryogenic gases. They protect against contact with cold objects and against liquefied gas splashes (liquid nitrogen). The hazards deriving from handling cryogenic liquefied gases are in direct relation to the extremely low temperatures of such substances. Exposing skin to extremely low temperatures can cause damage similar to burns. Although the gloves withstand contact with liquid nitrogen, contact with such a substance must be accidental: do not immerse the gloves in liquefied gas. Contact with liquefied gas hardens the materials used to make the gloves: in case of contact check after 20/30 seconds that the gloves have returned to their original softness and that there are no cracks/holes.

RISKS

The gloves are suitable against risks of:

- abrasion; small cut by blade/sheet metal;
- tearing; puncture/penetration;
- contact with cold and cryogenic gases. The gloves are NOT suitable against risks of:
- heat and fire;
- chemicals;
- all risks not mentioned in the Information Note herein.

ALLERGENS

- The manufacturer is not aware of the presence of any allergens. Kindly report any cases observed of hypersensitivity or allergic reaction.
- Upon contact with the skin on particularly sensitive people, any glove could cause allergic reactions not envisaged by the manufacturer. In such cases, we recommend you seek immediate medical advice.

IDENTIFICATION AND CHOICE OF SUITABLE GLOVES

The choice of suitable glove model must be made according to the specific requirements of the workplace, the type of risk and the related environmental conditions. The responsibility for identifying and choosing the appropriate/suitable pair of gloves (PPE) lies with the employer. Consequently, it is a good idea to check, before use, the suitability of the characteristics of this glove model for your specific requirements.

PRELIMINARY CHECKS AND USE: WARNINGS

Before use, perform a visual inspection of the gloves to make sure they are in perfect condition, clean and intact. Should the gloves not be intact (visible damage such as loose seams, breakages or smudges), they must be replaced.

MAINTENANCE INSTRUCTIONS AND EXPIRY DATE

New gloves are packed in series of 5 pairs, in a packet containing the information note. The gloves must be kept in their original packet, in a dry place away from sources of heat. Avoid contact with solvents which could alter the characteristics.

In particularly strenuous conditions of use or in settings with special situations, the gloves may be subject to sudden and unexpected deterioration not envisaged by the manufacturer. Consequently, it is not possible to establish an "expiry date".

CLEANING

NO WASHING IS PERMITTED to avoid the glove losing its safety characteristics. After use, the exterior of the gloves should be cleaned with a light detergent solution and left to dry in the open. Any contaminants must immediately be removed.

Explanation of **EN388** levels

- abc dx
- (a) Abrasion resistance (values from 1 to 4)
 - (b) Blade cut resistance (values from 1 to 5)
 - (c) Tear resistance (values from 1 to 4)
 - (d) Puncture resistance (values from 1 to 4)
 - (e) Cur resistance (ISO 13997:1999)

Explanation of **EN511** levels

- abc
- (a) Convective cold level (from 1 to 4)
 - (b) Contact cold level (from 1 to 4)
 - (c) Water penetration:
(0=not waterproof, 1=waterproof)

ALL MODELS UNDERGO SPECIFIC TEST FOR CONTACT WITH LIQUID NITROGEN (-195.82 °C) FOR ONE MINUTE AND SUBSEQUENT BENDING TEST.

Kora Srl – Cryokit Cryogenic Protection
MilanoFiori – Strada 6 – Edificio A – Scala 13
Tel. 02 48841819 – info@cryokit.net
www.cryokit.net

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Storage / Storage, Use, Cleaning, Maintenance, and Disinfection Instructions

The apron and gaiters, both new and used, must be examined in advance, before each use, in order to ensure that there are no conditions of wear, abrasions, lacerations or cuts, such as to compromise the degree of protection required. Any contaminants must be immediately removed. Any contamination, alteration or incorrect use of the device dangerously reduces its performance. Store in a cool dry place away from light.

The results of the tests on the prototypes determined a level / class of protection provided by the PPE

PPE class 3

Accessories that can be used together with PPE

Gloves, aprons and visors

Protection classes referring to the different risk levels and the corresponding limits of use.

The apron and leggings do not protect against:

- Chemical agents other than those indicated
- Ionizing radiations
- Biological agents
- Any risk other than that indicated in this information note

Duration and / or expiry of the PPE / component (month / year), if applicable

Each apron or gaiter is contained in a polythene envelope with this information note inside. To waste of the envelope refer to the provisions of the municipality to which it belongs.

Risks against which the PPE is designed

Processes that require protection against splashes of liquid nitrogen and cryogenic gases in general in clinical laboratories, the pharmaceutical industry, cryopreservation, bio-medical, biological banks, the food industry, handling of dry ice, etc. Possible leakage of liquid in the transfer operations (splash).

Name, address and identification number of the notified body in charge of assessing the PPE (release of Forms B, C2, D) *

ANCCP CERTIFICATION AGENCY srl - Via dello Struggino, 6
I-57121 Livorno (Italy) Tel: +39.0586.209006 Fax: +39.0586.278450
e-mail: anccp@anccp.it - n° Notified Body 0302
Harmonized standards / techniques used
EN388: 2016 / EN511: 2006

Website (where the EU Declaration of Conformity can be accessed)
www.cryokit.net

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EN USER INSTRUCTION

BEFORE USING THIS PPE, THE USER MUST:

- Make sure that it is suitable for the intended use (by thoroughly reading these notes) and make sure that the ocular / visor component is devoid of scratches and / or cracks or defects that could negatively affect resistance to impacts and / or correct viewing.
- In order to reduce risk, if necessary replace the damaged ocular.
- In case the PPE is supplied disassembled, the user must assemble each component according to any instructions reported either on the packaging or as contained on a separate instruction leaflet to ensure that performances indicated by markings are achieved;
- Position any removable and/or repositionable ocular/visor in the envisaged operating configuration so that the protection intended by the markings will be guaranteed;
- Assemble components and/or accessories supplied using any operator's manual provided with the single model if specified.
- Adapt any adjustment systems. PPE must not accidentally move during use. Wear the PPE before exposure to risk conditions against which it is meant to safeguard and use it only for those fields of use as envisaged by the markings;
- PPE must not be exchanged between several users but are specifically designed for personal protection. In case during use the PPE is damaged as a result of an impact, burning and/or exposure to corrosion, the PPE will have to be replaced. The only admissible fields of use are those resulting from markings applied to the PPE, the meaning of which is expressed hereafter in subsequent paragraphs. In order to understand product performances - Would you please refer also to the table expressing scale numbers as a function of shielding properties of the lenses/filters against optical radiation and the table hereafter reported describing the level of impacts sustained by the PPE.

Marking	Applicability	Definition	Striking element	Impact speed
S	All types of eye protectors	Increased robustness	Steel ball Ø 22mm Min. weight of 43g	approx. 5,1m/s

Marking	Applicability	Definition	Striking element	Impact speed
F	All types of eye protectors	Low-energy impact		45m/s
B	Goggles	Medium energy impact	Steel ball Ø 6mm weight 0.86g	120 m/s
A	Face-shields	High-energy impact		190 m/s

Any filter performance is described in the markings by the scale number as combination of code number and of shade number joined by a dash "-": for welding filters only shade number exists.

Standards by fields of use - [Lens Marking]	Code No.	Shade No.
EN169:2002 on Welding Filters	2	1,2 to 16
EN170:2002 on Welding Filters	2	1,2 to 5
EN171:2002 on Infrared Filters	4	1,2 to 10
EN172:1994+A1:2000+A2:2001 on Sunglare Filters for Industrial Use	6	1,1 to 4,1

C = after code number indicates oculars with higher colour identification performance (re EN170:2002) - Articles marked with scale numbers ranging from 2-1,2 to 2-5 may alter colour recognition when worn unless they are marked 2C. T = indicates that the PPE was tested for protection against high speed particles at extreme temperatures -5°C and +55°C (FT or BT or AT). If the device does not bear the letter T, it must only be used against high-speed particles at room temperature.

Description of Fields of Use expressed by markings applied

Marking Meaning

None Basic Use - Protection against mechanical hazard and hazards originated from exposure to non ionizing radiation - laser excluded

3 Protection against liquids - Liquid spray or drops

4 Protection against large dust particles - particles > 5µm

5 Protection against gas and fine dust particles - particles < 5µm

8 Protection against electric arc generated from short circuit

9 Protection against projections of molten metals and hot-solids

H PPE designed to fit to small-size head

R Increased reflectance of the ocular in the IR

K Resistance of the ocular to surface damage by fine particles

N Resistance of the ocular to fogging

v Replacement ocular

Marking of PPE in which frame and ocular form a single unit: complete ocular marking "-" of separation, EN166:2001 reference and pertinent symbols for fields of use and impact level.

Important - Compatibility of Markings: if oculars and frame do not bear the same symbols whether these are S, F or B or A or S, FT or BT or AT the entire equipment is downgraded to the lower class of impact resistance in respect of the field of use identified by markings.

Markings - Specifications: oculars classified in optical class 3 are not suitable or made for prolonged use. For conformity to the field of use identified by symbol 8, the PPE concerned must be equipped with a filter characterized by a Scale No. of either 2-1,2 or 3-1,2 and by a minimum thickness of 1,4mm. For specification on PPE assessment in respect of protection against electric arc on the basis of GS-ET-29 standard, refer to complementary user instruction. For conformity with symbol 9, frame and eyepiece must be marked 9 along with one of the symbols B or A. The eye protectors against high-speed particles worn above the optical eyepiece can transmit impacts and therefore generate hazards to the wearer. Filters showing the symbol are not suitable for road and driving use (limitation also applicable to filters with marking 5-4,1 and 6-4,1 as sunglare filters according to EN172:1994+A1:2000+A2:2001). For performances recorded in respect of this article, please refer to the marking applied to it and to the explanation of their meaning as described on to this leaflet. Risks against which this article is intended to provide protection against are expressed by markings applied to it.

The code of identification of this PPE is shown on the package. Article validated on the basis of EN166:2001.

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[] Standard number (EN166 or 166)

[] Symbols for pertinent field's of use

[] Symbol of increased robustness / resistance to high-speed particles / extreme temperatures (if pertinent)

[] Symbol indicating that the PPE is dimensionally designed for a small-size head. Symbol H (if pertinent)

[] Max scale number/s of the oculars - compatible with article's frame (if pertinent)

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[] Optical class (with the exception of covering shields) Symbol 1 or 2 or 3

[] Symbol of mechanical strength: Depending on the case S or F or B or A

[] Symbols representative of fields of use intended by the manufacturer

[] Symbol of resistance to surface damage by fine particles: K (if performance is guaranteed)

[] Symbol of resistance to fogging of oculars: N (if performance is guaranteed)

Specification according to EN175:1997 - This standard applies to PPE for eye and face protection during welding and related procedures. The conforming PPE is designed to include protection filters with or without protective/reinforced oculars, certified in conformity to EN166:2001 and to EN169:2002.

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EN175 Reference to EN175 standard

[] Symbols for the fields of use (if pertinent)

[] Weight in grams (when pertinent)

PPE must not be used in contexts and exposure conditions other than any expressed by markings or reported on the packs. Oculars / filters made of reinforced mineral glass can only be used associated with adequately reinforced ocular. Maintenance by the user is limited to cleaning procedure unless otherwise indicated on the package or in other complementary instruction leaflet. Any chat made to the PPE not performed by Univet as well as, any maintenance of components and/or replacement operation not performed in conformity with Univet instructions, shall result in the inva-

lidation of performances and of product type approval. Spare parts can be obtained from Univet s.r.l. Replacement of product structural parts, unless illustrated on the packs or on separate notes shall be made by Univet s.r.l., which reserves the right to decide as regard to complete replacement or article repair. Articles affected by wear, use not in compliance with specifications or tampered with shall be excluded from maintenance and/or replacement. Keep non-contaminated PPE in their original pack at a temperature between 5°C and 40°C in a dry place where there are no organic vapours or abrasive dusts in suspension. Store potentially contaminated PPE in a special container different from the supply pack and keep these separated from uncompromised devices until disposal or disinfection. To prevent any drop in product performances, do not expose to heat sources or UV radiations. The expiry date of the PPE depends on the context of use. If instructions for use, maintenance, cleaning and storage of PPE are followed, the article can be used for up to 3 years from when the pack is first opened. Process in any case to replace the PPE every 3 years and even earlier if evident signs of degradation appear. Clean PPE every time they have been used before returning them to the pack and check the integrity of each component. PPE can be cleaned with water and dried with a soft non-abrasive cloth; lenses with anti-fogging treatment must only be cleaned using a soft cloth.

For disinfection it has to be noted that different contaminating elements will require different type of approach and of disinfectants. As a general rule, Univet suggests to:

- Avoid to use any disinfectant containing: alcohol / ammoniaca / benzalconio cloride / chlorine / benzene / sodium hypochlorite/ solvents / acid and alkaline mixtures.
- Avoid any kind of abrasive paste - Use a mild detergent soap while taking care not to transfer it on the optical portion of the lenses. After deteration with soap, luke care to rinse the frame with cool water and to dry it with a soft non abrasive cloth.
- If the nature of contaminants will pose a risk to the safety and health even after completion of cleaning procedures, it is advisable to replace the entire PPE. Susceptible individuals may suffer from allergic reaction if some materials of which this article is made of will come into direct contact with their skin.

If this symbol is represented on the label applied to the package, latex is present in at least one component of this PPE.

Product made in Taiwan