## NOTA INFORMATIVA cod. 393061 - mod. Nitrile 55S

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## INFORMATIVE NOTE cod. 393061 - mod. Nitrile 55S

Read this briefing note through carefully prior to use and prior to any maintenance work. The purpose of the information contained herein is to assist and direct the user in choosing and using the equipment.

Neither the manufacturer nor the distributor can be held liable for an incorrect use of the equipment.

This briefing note must be kept for the duration of the equipment.

Glove model: **cod. 393061 - mod. Nitrile 55S.**

## CE2777

The CE 2777 mark means that these gloves are personal protective equipment conforming to the essential health and safety requirements set forth in the Regulation (UE) 2016/425 and have been certified by the notified body: **SATRA Technology Europe Ltd, Bracetown Business Park, Clonee, Dublin, D15 YN2P, Ireland n° 2777**.

The manufacturer is subject to the control of CE quality guarantee of the production (Module D) by the notified body: **SATRA Technology Europe Ltd, Bracetown Business Park, Clonee, Dublin, D15 YN2P, Ireland n° 2777**.

OUR SIZE REFERENCES	SIZES CONTEMPLATED EN ISO 21420:2020
S	7
M	8
L	9
XL	10

The dexterity level of the fingers is equivalent to 5 out of 5.

The gloves meet the safety characteristics only if the size is right and worn properly.

#### RANGE OF USE

The glove is designed to protect the user against contact with chemical products and micro-organisms (intended as fungus, bacteria and viruses).

The glove can be used also as a CLASS I MEDICAL DEVICE, pursuant to the legal provisions of Regulation 2017/745/EU, for use in the professional sector and within the structures of the National Health Service, provided the warnings on use and care indicated on the glove packaging are observed.

Gloves certified according to the EN 455 standard.

#### PICTOGRAMS AND PROTECTION LEVELS

 <b>EN ISO 374-1:2016 +A1:2018 TYPE B</b>	 <b>EN ISO 374-5:2016</b>
<b>KPT</b>	<b>VIRUS</b>

#### PROTECTION AGAINST CHEMICALS

**PROTECTION AGAINST MICRO-ORGANISMS**

Resistance to air and water penetration using the method referred to in EN ISO 374-2:2019 standards

CHEMICAL RESISTANCE TABLE according to EN ISO 374-1:2016+A1:2018


Chemical	Permeation index (0 - 6)	Time of Passage (min)	Degradation index (EN ISO 374-4:2019)
K Sodium hydroxide 40%	6	480	-68.1%
P Hydrogen peroxide 30%	2	30	30.5%
T Formaldehyde 37%	5	240	9.5%

Time of passage.

The time a chemical takes to pass through the glove material.

The values may differ when the chemical substance is used in a compound.

Do not use the gloves for handling the above-mentioned chemical products for longer than the times stated above.

 <b>Class I medical device</b> (rule 5, Chapter III, Annex VIII Regulation 2017/745/EU) disposable, non-sterile suitable for general use (exploration and hygiene).			
<b>COMPLIANT WITH: REGULATION (EU) N. 10/2011 - REGULATION (EC) N. 1935/2004 - DIRECTIVES 85/572/EEC AND SUBSEQUENT AMENDMENTS.</b>			
Ethanol 10%	Acetic acid 3%	Vegetable oil	Polyphenyl oxide (MPPO)
Yes	No	Yes	No
30 min. / 40°C	-	30 min. / 40°C	-

Since the performance levels are based on lab test results, they may differ from the actual conditions at the workplace as temperature, abrasion and degradation. Consequently it is advisable to conduct a preliminary test of the gloves to see if they are suitable for their intended use.

The protective characteristics refer to new equipment, in good condition and that has never been subject to cleaning treatments.

Motions, tear, frictions and degradations caused by the contact with a chemical substance could reduce significantly the effective usage of the frame in the case of corrosive chemicals, degradation may be the most important factor to consider when choosing gloves.

#### WARNINGS

Prior to use, make sure the gloves are in good condition by which we mean no cuts, holes or cracks, etc. If the gloves are not in perfect condition, replace them immediately. The gloves must be used solely for the risks contemplated in this briefing note. Avoid using PPE near moving parts in which it could get entangled.

**Hand Hygiene**: When an indication of hand hygiene precedes a contact that also requires glove usage, hand rubbing or hand washing should be performed before donning glove.

**Glove donning and doffing**: When donning the glove, please ensure that both glove and hand are clean, glove size is right, and it is properly fit on the contours and crotches of fingers, in order to avoid the transmission of infections. When doffing the glove, avoid letting the external surface of the gloves come into contact with the skin, as the surface could be contaminated with blood or other body fluids.

After doffing, powder can be removed by wiping gloves thoroughly with a sterile wet sponge, sterile wet towel or other effective methods.

Avoid breakage as it may result in splashes of contaminants into the eyes, mouth, skin or above other nearby subjects.

**If used as a NON STERILE medical device, its use must be transient if used invasively, in relation to the natural orifices of the body. IT MUST NOT be used in contact with open wounds, or for surgery or inside unnatural orifices of the body or wherever the use of a STERILE device is requested.**

#### STORING

The gloves are packed in cardboard boxes.

Keep the gloves in their original packaging in a clean, dry place, away from sources of heat and direct sunlight. The user must always inspect the gloves to ensure they are sound prior to using them. The glove must be used, provided the packaging is perfectly intact and kept in the preservation conditions stated on the packaging itself, **within the indicated expiry date and no later**. Since these are disposable gloves, they can be used just once; they must be replaced at the end of their use - they cannot be used a second time because they could be contaminated and there may be risk of infection or spread of pathogens.

#### NOTES

In case of manufacturing defects, the device will be replaced.

The used gloves might be contaminated by chemical or biological products, must be disposed in compliance with the local regulations (dump, incinerator). Latex, powder and DEHP/DOP free. Please report any cases found of hypersensitivity or allergic reactions. The EU conformity declaration is available at the following address: [www.nerispa.com](http://www.nerispa.com)

Please also report to the manufacturer and to the authorities any serious incident related to the medical device.

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